



# GSTPur Glutathione kit

## Index

|                                |   |
|--------------------------------|---|
| 1. Product Description.....    | 1 |
| 2. Purification Procedure..... | 1 |
| 3. Cleaning-in-Place.....      | 2 |
| 4. Troubleshooting.....        | 2 |
| 5. Related Products.....       | 3 |

## 1. Product Description

Glutathione Beads is an affinity chromatography medium designed for one-step purification of recombinant glutathione S-transferase (GST) fusion proteins and other glutathione binding proteins expressed in *E. coli*, insect cells and mammalian cells. The recombinant GST fusion proteins can be purified directly from pre-treated cell lysate using Glutathione Beads. It is the excellent choice for high performance purifications. **GSTPur Glutathione kit** contains one GSTPur Glutathione Column, Lysis Buffer and Elution Buffer. It is fast, simple and easy operation.

Table 1. Characteristics of Glutathione Beads

| Item                    | Description                               |
|-------------------------|-------------------------------------------|
| Matrix Spherical        | 4% agarose                                |
| Static Binding Capacity | >20mg GST-tagged protein(40kDa)/ml medium |
| Particle size           | 45-165um                                  |
| Maximum Pressure        | 0.1MPa, 1bar                              |
| Storage Solution        | 1X PBS containing 20% ethanol             |
| Storage Temperature     | 2-8℃                                      |

Table 2. The component of GSTPur Glutathione Kit

|                           | Cat.No   | Size  | Amount |
|---------------------------|----------|-------|--------|
| GSTPur Glutathione Column | SA008K-1 | 3ml   | 1      |
| Lysis Buffer              | SA008K-2 | 100ml | 2      |
| Elution Buffer            | SA008K-3 | 100ml | 1      |

## 2. Purification Procedure

### 2.1 Sample Preparation

It is recommended to filter the sample solution by passing them through a 0.22µm or 0.45 µm filter before use.

### 2.2 Column Purification

**GSTPur Glutathione kit** contains GSTPur Glutathione Column, Lysis Buffer and Elution Buffer. The purification process takes about 30min. It mainly depends on the sample volume and solution viscosity.

- 1) Fix GSTPur Glutathione Column and remove the bottom plug and top plug. Allow storage buffer to drain from resin by gravity flow.
- 2) Equilibrate column with 5ml Lysis Buffer. Allow buffer to drain from the column. Repeat this step two times.
- 3) Add the prepared protein extract to the resin. Collect the flow-through. If desired, re-apply the flow-through once to maximize binding.
- 4) Wash resin with 5ml Lysis Buffer and collect the flow-through. Repeat this step five times or repeat this step using a new collection tube until the absorbance of the flow-through fraction at 280nm approaches baseline.
- 5) Elute GST-tagged proteins from the resin with 15-30ml Elution Buffer. Collect the fraction every 5ml in a separate tube.
- 6) Equilibrate the column with 5ml Lysis Buffer and 5ml deionized water alternately. Repeat 2 times. Then Equilibrate the column with 5ml 20% ethanol. Finally store the beads with 1XPBS containing 20% ethanol at 2-8℃.

### 2.4 Analysis

Identify the fractions containing the protein. Use UV absorbance, SDS-PAGE, or western blot.






Fig.2. Purification process of GSTPur Glutathione Column

Step1: Equilibrate; Step2: binding; Step3: Wash; Step4: Elution

### 3. Cleaning-in-Place

GSTPur Glutathione Column can be reused to purify the same protein three times without regeneration. If the target GST-fusion protein is different, however, the Glutathione resin must be regenerated using the following protocol:

#### Remove the precipitation or denatured protein

Wash the column with 2 column volumes 6M guanidine hydrochloride solution. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

#### Remove the hydrophobically bound protein

Wash the column with 3-4 column volumes 70% ethanol or 2 column volumes 1% Triton X-100. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

### 4. Troubleshooting

| Problem                                                              | Probable Cause                                                                | Solution                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------------------------------------------------------|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The yield of the purified GST fusion protein is low or undetectable. | The fusion protein forms inclusion body.                                      | Grow bacteria at lower temperature (20-30°C), or reduce final concentration of IPTG to 0.1 mM for protein induction, or reduce the induction time.<br>Properly dissolve and refold the inclusion body prior to the purification.                                                                                                                                              |
|                                                                      | The fusion protein does not bind to Glutathione Beads efficiently.            | Use batch method for purification. Incubate clear solution (the sonicate, etc) containing GST-fusion protein with Glutathione Beads for 2 hours or longer (such as overnight) and then load the mixture onto the column.                                                                                                                                                      |
|                                                                      | The fusion protein does not contain active GST.                               | Use mild sonication condition or other lysis method, such as lysozyme so that GST is not denatured.                                                                                                                                                                                                                                                                           |
|                                                                      | The fusion protein is degraded by protease.                                   | Add appropriate protease inhibitors such as PMSF in the lysis solution and wash solution.                                                                                                                                                                                                                                                                                     |
|                                                                      | The fusion protein is not efficiently eluted from Glutathione Beads.          | Increase elution time or increase the concentration of glutathione to 15 mM or higher in the elution buffer.<br>Adjust the pH of the elution buffer to 8.0-9.0 without increasing the glutathione concentration.<br>Add Triton X-100 (0.1%, final concentration) or Noctylglucoside (2%, final concentration) or NaCl (0.1-0.2 M, final concentration) to the elution buffer. |
| Multiple bands observed in the eluted protein                        | The fusion protein is degraded by protease.                                   | Add appropriate protease inhibitors (or inhibitor cocktails) such as PMSF in the lysis solution and wash solution.                                                                                                                                                                                                                                                            |
|                                                                      | Some host proteins,such as chaperonins, may interact with the fusion protein. | Add DTT (5 mM, final concentration) in the wash buffer. Incubate the recombinant protein solution in chaperonin buffer (2 mM ATP, 10 mM MgSO <sub>4</sub> , 50 mM Tris-HCl) at 37°C for 10 min prior to the purification.                                                                                                                                                     |
|                                                                      | Over-sonication will cause some protein to bind to the fusion protein.        | Use milder sonication condition or another lysis method.                                                                                                                                                                                                                                                                                                                      |
|                                                                      | Some protein will bind to the fusion protein or beads non-specifically.       | Optimizing the wash conditions. Detergents such as 1% Triton X-100, 1% Tween-20, 0.03% SDS, or 0.1% NP-40 may be used to reduce non-specific binding. Salt concentration in the wash solution can also be optimized to reduce non-specific binding.                                                                                                                           |





## 5. Related Products

| Product                | Cat. No. | Size          |
|------------------------|----------|---------------|
| Glutathione Beads      | SA008005 | 5 ml          |
|                        | SA008025 | 25 ml         |
|                        | SA008100 | 100 ml        |
|                        | SA008500 | 500 ml        |
|                        | SA00801L | 1 L           |
|                        | SA00810L | 10 L          |
| GSTPur Glutathione Kit | SA008K03 | 3 次           |
|                        | SA008K05 | 5 次           |
| Glutathione Beads 4FF  | SA010005 | 5 ml          |
|                        | SA010025 | 25 ml         |
|                        | SA010100 | 100 ml        |
|                        | SA010500 | 500 ml        |
|                        | SA01001L | 1 L           |
|                        | SA01010L | 10 L          |
| GSTCap 4FF             | SA010C11 | 1X1 ml        |
|                        | SA010C51 | 5x1 ml        |
|                        | SA010C15 | 1X5 ml        |
|                        | SA010C55 | 5X5 ml        |
|                        | SA010CS  | 3X1 ml+1X5 ml |

