



# HisCap IDA 6FF

## Index

|                                |   |
|--------------------------------|---|
| 1. Product Description.....    | 1 |
| 2. Purification Procedure..... | 1 |
| 3. Cleaning-in-Place.....      | 3 |
| 4. Regeneration.....           | 3 |
| 5. Troubleshooting.....        | 3 |
| 6. Related Products.....       | 4 |

## 1. Product Description

**Ni IDA Beads 6FF** can be used to purify 6xHis-tagged proteins expressed in series of expression systems, such as *E.coli*., yeast, insect cells and mammalian cells. **Ni IDA Beads 6FF** consists of 90µm beads of highly cross-linked 6% agarose, to which iminodiacetic acid (IDA) has been coupled. The chelating group has then been charged with nickel ions ( $\text{Ni}^{2+}$ ). This form is relatively stable plane structure, which can be destroyed by other competitive small molecule. The structure of Ni-IDA has more sites to bind histidine-tagged proteins. **Ni IDA Beads 6FF** has the advantages of higher binding capacity and cost-effective. Its high flow properties make it excellent for scale-up. See table 1.

**HisCap IDA 6FF** is one of a range of prepacked, ready-to-use columns for affinity chromatography. It is packed with 1ml and 5ml of **Ni IDA Beads 6FF**. Five different packing sizes are available. **HisCap IDA 6FF** has the standard interface and can be adapted to all kinds of chromatography system, such as ÄKTA. It is fast, simple and easy operation.

Table 1. Characteristics of Ni IDA Beads 6FF

| Item                    | Description                           |
|-------------------------|---------------------------------------|
| Matrix Spherical        | Highly cross-linked 6% agarose        |
| Static Binding Capacity | >40 mg 6xHis-tagged protein/ml medium |
| Particle size           | 45-165µm                              |
| Maximum Pressure        | 0.3MPa, 3bar                          |
| Storage Solution        | 1xPBS containing 20% ethanol          |
| Storage Temperature     | 4-30℃                                 |

## 2. Purification Procedure

### 2.1 Buffer Preparation

The basic principle of the following recommended buffer and other buffer is low concentration of imidazole in Lysis and wash buffer and high in elution buffer. Water and chemicals used for buffer preparation should be of high purity. It is recommended filtering the buffers by passing them through a 0.22µm or 0.45 µm filter before use. The recommended buffer is shown in table 2.

Table 2. Recommended buffer of his-tigged protein purification under native conditions

| Name           | Volume | Ingredient                                                                                                                                                                                                            |
|----------------|--------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Lysis Buffer   | 1L     | 50 mM $\text{NaH}_2\text{PO}_4$ (7.80 g $\text{NaH}_2\text{PO}_4 \cdot 2\text{H}_2\text{O}$ )<br>300 mM NaCl (17.54 g NaCl)<br>10 mM imidazole (0.68 g imidazole)<br>Adjust the buffer pH to 8.0 with NaOH solution.  |
| Wash Buffer    | 1L     | 50 mM $\text{NaH}_2\text{PO}_4$ (7.80 g $\text{NaH}_2\text{PO}_4 \cdot 2\text{H}_2\text{O}$ )<br>300 mM NaCl (17.54 g NaCl)<br>20 mM imidazole (1.36 g imidazole)<br>Adjust the buffer pH to 8.0 with NaOH solution.  |
| Elution Buffer | 1L     | 50 mM $\text{NaH}_2\text{PO}_4$ (7.80 g $\text{NaH}_2\text{PO}_4 \cdot 2\text{H}_2\text{O}$ )<br>300 mM NaCl (17.54 g NaCl)<br>250 mM imidazole (17.0 g imidazole)<br>Adjust the buffer pH to 8.0 with NaOH solution. |

### 2.2 Sample Preparation

#### 2.2.1 Protein expressed in *E.coli* or yeast

1) Single colonies were cultured in LB medium. According to the instruction, adding the inducers for a period of time.





- 2) Harvest cells from an appropriate volume of bacterial culture by centrifugation at 7,000 rpm(7,500×g) for 10-15 min at 4℃. Discard supernatant and determine weight of pellet. Resuspend pellet in 1:10 ration(w/v) in lysis buffer and add lysozyme (0.2-0.4mg/ml cell paste ,if the host cell containing pLysS or pLysE ,it can be without lysozyme ) and PMSF (1mM/ml cell paste) .
- 3) If high concentration of cell suspension ,it is consider to add 10 µg/ml RNase A 和 5 µg/ml DNase I .Sonicate the cell suspension/lysate on ice.
- 4) Centrifuge the homogenized lysate at 10,000 rpm(15,000×g) for 20 min at 4℃ to clarify sample. Save supernatant.

## 2.2.2 Protein expressed in yeast, insect or mammalian cells

- 1) Harvest the cells from an appropriate volume of culture by centrifugation at 5,000 rpm(3,800×g) for 10-15 min at 4℃. Save supernatant. If the supernatant without EDTA, histidine and reductant, it can be purified directly, otherwise it need dialysis to 1XPBS under 4℃ .
- 2) For a large volume of supernatant, it need precipitation by adding ammonium sulfate and dialysis to 1XPBS under 4℃ .

## 2.3 Sample Purification

- 1) Fill the syringe or pump tubing with distilled water. Remove the stopper and connect the column to the syringe (with the provided connector), or pump tubing, "drop to drop" to avoid introducing air into the column. Remove the snap-off end at the column outlet.
- 2) Wash the column with 3-5 column volumes of distilled water.
- 3) Equilibrate the column with at least 5 column volumes Lysis Buffer.
- 4) Apply the pre-treated sample, using a Loop fitted to the connector or by pumping it onto the column.
- 5) Wash with Wash Buffer until the absorbance reaches the baseline or no material appears in the effluent(Generally at least 10-15 column volumes).
- 6) Elute with elution buffer using a stepwise or linear gradient. For one-step elution, 5 column volumes are usually enough. Other volumes may be required if the interaction is difficult to break. Linear gradient elution can be used to separate proteins of different binding strengths with a small gradient, such as 20 column volumes or more.

## 2.4 Analysis

Identify the fractions containing the His-tagged protein. Use UV absorbance, SDS-PAGE, or western blot.

## 3. Cleaning-in-Place

A column used to purify protein from cell extract usually contains some soluble substances and cell debris that are nonspecifically absorbed onto the matrix. Cleaning-in-Place eliminates material not removed by regeneration and prevents progressive buildup of contaminants. If the column is to be reused ,these contaminants should be cleaned from the column, as they were not completely removed during the sample clarification steps.

### ● Remove the strong hydrophobic binding protein, lipoprotein and lipid

Wash the column with 5-10 column volumes 30% isopropanol contacting for 15-20 min. Or you can choose the 2 column volumes acidic or alkaline solution containing detergents, for example, 0.1 M acetic acid solution contains 0.1-0.5% non-ionic detergent, contacting for 1-2 hours .Finally wash the column with 10 column volumes distilled water.

### ● Remove the proteins combined with ion interacting

Wash the column with 1.5M NaCl solution contacting for 10-15 min.Finally wash the column with 10 column volumes distilled water.

## 4. Regenerating the medium

In general, **Ni IDA Beads 6FF** may be reused before it becomes necessary to recharge them with metal ions. When the back pressure is too highor the capacity significantly lower, it need to strip the metal ions and recharge the **Ni IDA Beads 6FF** as the following procedure. Wash the column with one of the following solutions.

- 1) 5 column volumes of distilled water;
- 2) 5 column volumes of 100 mM EDTA (pH 8.0);
- 3) 10 column volumes of distilled water;
- 4) 5 column volumes of 0.5 M NaOH , contacting for 10-15 min;
- 5) 10 column volumes of distilled water;
- 6) 3-5 column volumes of 100 mM NiSO<sub>4</sub> ;
- 7) 10 column volumes of distilled water;

After regeneration, the medium can be used immediately, otherwise, it need to be suspended in an equal volume of 1X PBS containing 20% ethanol at 4℃.





## 5. Troubleshooting

| Problem                                  | Probable Cause                                                             | Solution                                                                                                                                                              |
|------------------------------------------|----------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Back pressure exceeds 3 bar              | Column is clogged                                                          | Cleaning in place(part 3).<br>Increase the centrifugation speed or filtering the sample.                                                                              |
|                                          | Sample is too viscous                                                      | Increase sonication or add DNase I (5 µg/ml with 1mM Mg2+. Incubate on ice for 15min.                                                                                 |
|                                          | Buffer is too viscous                                                      | Dilute sample by adding more homogenization buffer.                                                                                                                   |
| No protein is eluted                     | Expression of target protein in extract is very low                        | Check expression level of protein by estimating the amount in the extract, flow through, elute fraction and pellet upon centrifugation. Or apply large sample volume. |
|                                          | Target protein is found in the flow through                                | Reduce imidazole concentration in lysis buffer sample and wash buffer. Or increase buffer pH.                                                                         |
|                                          | Elution conditions are too mild.                                           | Increase imidazole concentration in Elution buffer. Or decrease buffer pH.                                                                                            |
|                                          |                                                                            | Strip nickel ion by using 10-100mM EDTA solution, at the same time you can obtain target protein.                                                                     |
|                                          | Protein degradation or purification cause the histidine tag to be removed. | Operate at 4°C. Add protease inhibitors.<br>Make a new construct with his-tag attached to other terminus.                                                             |
| His-tagged protein is not pure           | Wash is not enough                                                         | Increase the volume of Wash Buffer.                                                                                                                                   |
|                                          | Association between the his tagged protein and protein contaminant.        | Optimize the wash condition by adjusting the pH and imidazole concentration.                                                                                          |
|                                          |                                                                            | Add an additional chromatography step, that is ion exchange, hydrophobic interaction or size exclusion.                                                               |
| The color of medium becomes shallow.     | The nickel ions was stripped.                                              | Chelate nickel ions according to the part 4.                                                                                                                          |
| The medium are brown.                    | The buffer contains DTT.                                                   | Use Ni NTA Beads or Ni NTA Beads FF when the reductant concentration bellow 2 mM.                                                                                     |
| Protein precipitates during purification | Temperature is too low                                                     | Perform the purification at room time.                                                                                                                                |
|                                          | Aggregate formation                                                        | Add solubilization agents to the samples and buffers, for example 0.1% Triton X-100 , Tween-20 and ≤20% glycerol to maintain protein solubility.                      |

## 6. Related Products

| Product          | Cat. No. | Size          |
|------------------|----------|---------------|
| Ni IDA Beads     | SA003005 | 5 ml          |
|                  | SA003025 | 25 ml         |
|                  | SA003100 | 100 ml        |
|                  | SA003500 | 500 ml        |
|                  | SA00301L | 1 L           |
|                  | SA00310L | 10 L          |
| Ni IDA Beads 6FF | SA052005 | 5 ml          |
|                  | SA052025 | 25 ml         |
|                  | SA052100 | 100 ml        |
|                  | SA052500 | 500 ml        |
|                  | SA05201L | 1 L           |
|                  | SA05210L | 10 L          |
| HisCap IDA 6FF   | SA052C11 | 1×1 ml        |
|                  | SA052C51 | 5×1 ml        |
|                  | SA052C15 | 1×5 ml        |
|                  | SA052C55 | 5×5 ml        |
|                  | SA052CS  | 3×1 ml+1×5 ml |
|                  |          |               |

