



# PreCap Streptavidin ST

## Index

|                                 |   |
|---------------------------------|---|
| 1. Product Description.....     | 1 |
| 2. Purification Procedure ..... | 2 |
| 3. Regeneration.....            | 2 |
| 4. Related Products .....       | 3 |

## 1. Product Description

**Streptavidin ST Beads 4FF** is a chromatography medium for one-step purifying Strep-tag II fusion proteins from various of expression system. The Strep-Tag II peptide is an eight amino acid fusion tag(Trp-Ser-His-Pro-Gln- Phe-Glu-Lys), which has negligible effects on recombinant proteins. The ligand immobilized on high-crossed linked 4% agarose is a specially recombinant protein. Purification under physiological conditions and mild elution preserves the activity of the target protein.

**PreCap Streptavidin ST** is a prepacked ready to use column. **PreCap Streptavidin ST** 1ml and 5ml columns are packed with 1ml and 5ml of **Streptavidin ST Beads 4FF**. **PreCap Streptavidin ST** can be adapted to all kinds of chromatography system, such as ÄKTA. It is easy to operate.

Table 1. Characteristics of **Streptavidin ST Beads 4FF**

| Item                  | Description                          |
|-----------------------|--------------------------------------|
| Matrix                | Highly cross-linked 4% agarose beads |
| Ligand                | Streptavidin ST                      |
| Capacity (/ml medium) | 6 mg Strep-tag II fusion proteins    |
| Particle Size (μm)    | 45-165                               |
| Maxi Pressure         | 0.3 MPa, 3 bar                       |
| pH                    | 3-10                                 |
| Storage Buffer        | 1XPBS containing 20% ethanol         |
| Storage Temperature   | 2°C - 8°C                            |

Table 2. Chemical compatibilities for **Streptavidin ST Beads 4FF**

| Reagent                                 | Concentration |
|-----------------------------------------|---------------|
| <b>Reduction Agents</b>                 |               |
| DTT                                     | 50 mM         |
| β-mercaptoethanol                       | 50 mM         |
| <b>Non-Ionic Detergents</b>             |               |
| C8E4 Octyltetraoxyethylene              | Max.0.88 %    |
| C10E5; Decylpentaoxyethylene            | 0.12 %        |
| C10 E6                                  | 0.03 %        |
| C10E8                                   | 0.005 %       |
| C12E9; Dodecyl nonaoxyethylene (Thesit) | 0.023 %       |
| DM; Decyl-β-D-maltoside                 | 0.35 %        |
| LM; N-dodecyl β-D-maltoside             | 0.007 %       |
| NG; N-nonyl-β-D-glucopyranoside         | 0.2 %         |
| OG; N-octyl-β-D-glucopyranoside         | 2.34 %        |
| TX; Triton X-100                        | 2 %           |
| Tween-20                                | 2 %           |
| <b>Ionic Detergents</b>                 |               |
| N-lauryl-sarcosine                      | 2 %           |
| 8-HESO;N-octyl-2-hydroxy-ethylsulfoxide | 1.32 %        |
| SDS; Sodium-N-dodecyl sulfate           | 0.1 %         |





| Zwitter-Ionic Detergents                                         |            |
|------------------------------------------------------------------|------------|
| CHAPS                                                            | 0.1 %      |
| DDAO; N-decyl-N,N-dimethylamine-N-oxide                          | 0.034 %    |
| LDAO; N-dodecyl-N,N-dimethylamine-N-oxide                        | 0.13 %     |
| Other reagent                                                    |            |
| Ammonium sulfate (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> | 2 M        |
| CaCl <sub>2</sub>                                                | Max.1 M    |
| Ethanol                                                          | 10%        |
| EDTA                                                             | 50 mM      |
| Guanidine                                                        | Max.1 M    |
| Glycerol                                                         | Max.25 %   |
| Imidazole                                                        | Max.250 mM |
| MgCl <sub>2</sub>                                                | 1 M        |
| NaCl                                                             | 5 M        |
| Urea                                                             | Max.1 M    |

## 2. Purification Procedure

### 2.1 Buffer Preparation

Water and chemicals used for buffer preparation should be high purity. It is recommended to filter the buffers by passing them through a 0.22 µm or 0.45 µm filter before use.

**Binding/Wash buffer:** 100 mM Tris-HCl, 150 mM NaCl, 1mM EDTA, pH 8.0 or PBS

**Elution Buffer:** 2.5mM desthiobiotin in binding buffer

**Regeneration Buffer:** 1mM HABA in binding buffer

### 2.2 Sample preparation

The sample should be adjusted to the composition of the binding buffer. This can be done either by diluting the sample with binding buffer or by buffer exchange. It is recommended to filter the sample solution by passing them through a 0.22 µm or 0.45 µm filter before use.

### 2.3 Sample Purification

- 1) Fill the syringe or pump tubing with binding buffer. Remove the stopper and connect the column to the syringe (with the provided connector), or pump tubing, "drop to drop" to avoid introducing air into the column. Remove the snap-off end at the column outlet.
- 2) Wash the column with 10 column volumes of binding buffer.
- 3) Apply the sample, using a syringe fitted to the connector or by pumping it onto the column.
- 4) Wash with 5 to 10 column volumes of wash buffer or until no material appears in the effluent.
- 5) Elute with 5 column volumes of elution buffer. Other volumes may be required if the interaction is difficult to break.

### 2.4 Analysis

Identify the fractions containing the target protein. Use UV absorbance, SDS-PAGE, or western blot.

## 3. Regeneration

**Regeneration:** Wash the column with 3CV distilled water followed by 15CV regeneration buffer and 30CV binding buffer. The displacement is detected by the change in color of the edium in the column from yellow to red. This color change is due to the accumulation of HABA/Streptavidin ST complexes. The HABA is washed away with the binding buffer.

**Equilibration:** Before next use, balance with 5 times column volume of equalizer.





#### 4. Related Products

| Product                   | Cat. No  | Size          |
|---------------------------|----------|---------------|
| Streptavidin ST Beads 4FF | SA053005 | 5 ml          |
|                           | SA053025 | 25 ml         |
|                           | SA053100 | 100 ml        |
|                           | SA053500 | 500 ml        |
|                           | SA05301L | 1 L           |
|                           | SA05310L | 10 L          |
| PreCap Streptavidin ST    | SA053C11 | 1X1 ml        |
|                           | SA053C51 | 5x1 ml        |
|                           | SA053C15 | 1X5 ml        |
|                           | SA053C55 | 5X5 ml        |
|                           | SA053CS  | 3X1 ml+1X5 ml |

