



# HiPur A 4FF

## Index

|                                |   |
|--------------------------------|---|
| 1. Product Description.....    | 1 |
| 2. Purification Procedure..... | 2 |
| 3. Maintenance.....            | 2 |
| 5. Troubleshooting.....        | 3 |
| 6. Related Products.....       | 3 |

## 1. Product Description

Smart-Lifesciences **rProtein A Beads 4FF** is an affinity chromatography medium designed for the purification of monoclonal antibody, polyclonal antibody and Fc tag fusion protein. The recombinant protein A is coupled to highly cross-linked 4% sepharose gel. This coupling method can improve the binding capacity for immunoglobulin. The static binding capacity of **rProtein A Beads 4FF** is greater than 40 mg human IgG/ml settled beads. The dynamic binding capacity will vary depending on several factors such as target antibody, flow rate etc. The characteristics of **rProtein A Beads 4FF** are summarized in Table 1 .

Protein A ,a bacterial cell wall protein isolated from *Staphylococcus aureus*, combines mammalian IgGs mainly through Fc regions. Recombinant protein A contains five high affinity IgG binding domains with other non-essential domains removed to reduce nonspecific binding.

**HiPur A 4FF** is a prepacked ready to use column. It is packed with 20ml rProtein A Beads 4FF .The arrow on the column label indicates the recommended flow direction. HiPur A 4FF can be adapted to all kinds of chromatography system, such as ÄKTA. It is easy to operate.

Table 1. Characteristics of **HiPur A 4FF**

| Item                    | Description                          |
|-------------------------|--------------------------------------|
| Size                    | 20 ml                                |
| Column size             | 1.6×10 cm                            |
| Matrix Spherical        | Highly cross-linked 4% agarose beads |
| Ligand                  | recombinant protein A                |
| Static Binding Capacity | >40mg Human IgG/ml medium            |
| Particle size           | 45-165um                             |
| Maximum Pressure        | 0.3MPa, 3bar                         |
| pH                      | 3-10                                 |
| Storage Solution        | 1×PBS containing 20% ethanol         |
| Storage Temperature     | 2℃-8℃                                |

Table 2. Relative binding strengths of antibodies from various species to protein G and protein A as measured in a competitive ELISA test.

| Species        | Subclass | Protein A binding | Protein G binding |
|----------------|----------|-------------------|-------------------|
| Human          | IgA      | variable          | —                 |
|                | IgD      | —                 | —                 |
|                | IgE      | —                 | —                 |
|                | IgG1     | ++++              | ++++              |
|                | IgG2     | ++++              | ++++              |
|                | IgG3     | —                 | ++++              |
|                | IgG4     | ++++              | ++++              |
|                | IgGM     | variable          | —                 |
| Avian egg yolk | IgY      | —                 | —                 |
| Cow            |          | ++                | ++++              |
| Dog            |          | ++                | +                 |
| Goat           |          | —                 | ++                |
| Guinea pig     | IgG1     | ++++              | ++                |
|                | IgG2     | ++++              | ++                |
| Hamster        |          | +                 | ++                |





|                |                |          |      |
|----------------|----------------|----------|------|
| Horse          |                | ++       | ++++ |
| Koala          |                | —        | +    |
| Llama          |                | —        | +    |
| Monkey(rhesus) |                | ++++     | ++++ |
| Mouse          | IgG1           | +        | ++++ |
|                | IgG2a          | ++++     | ++++ |
|                | IgG2b          | +++      | +++  |
|                | IgG3           | ++       | +++  |
|                | IgM            | variable | —    |
| Pig            |                | +++      | +++  |
| Rabbit         | no distinction | ++++     | +++  |
| Rat            | IgG1           | —        | +    |
|                | IgG2a          | —        | ++++ |
|                | IgG2b          | —        | ++   |
|                | IgG3           | +        | ++   |
| Sheep          |                | +/-      | ++   |

++++=strong binding; +++ medium binding; —=weak binding or no binding

## 2. Purification Procedure

### 2.1 Buffer Preparation

Water and chemicals used for buffer preparation should be of high purity. It is recommended to filter the buffers by passing them through a 0.22 or 0.45 µm filter before use.

**Binding/Wash Buffer:** 0.15 M NaCl, 20 mM Na<sub>2</sub>HPO<sub>4</sub>, pH 7.0

**Elution Buffer:** 0.1 M glycine, pH 3.0

**Neutralization Buffer:** 1 M Tris-HCl, pH 8.5

### 2.2 Sample Preparation

To insure that proper ionic strength and pH are maintained for optimal binding, it is necessary to dilute serum samples, ascite fluid or cell culture supernatant at least 1:1 with Binding/Wash Buffer. Alternatively, the sample may be dialyzed overnight with Binding/Wash Buffer.

### 2.4 Sample Purification

**HiPur A 4FF** is a prepacked, ready to use column for purification. The prepacked column provides fast, simple and easy separations in a convenient format.

- 1) Fill the pump tubing with binding buffer. Connect the column to purification system, “drop to drop” to avoid introducing air into the column.
- 2) Wash the column with 10 column volumes of binding buffer .
- 3) Apply the sample, using a syringe fitted to the connector or by pumping it onto the column.
- 4) Wash with 5 to 10 column volumes of binding buffer or until no material appears in the effluent.
- 5) Elute with 5 column volumes of elution buffer. Other volumes may be required if the interaction is difficult to break.
- 6) Add 10µl Neutralization Buffer to each 100 µl of eluate to neutralize the pH.

### 2.5 Analysis

Identify the fractions using UV absorbance, SDS-PAGE, or western blot.

## 3. Maintenance

### 3.1 Regeneration

Regenerate the column by washing the resin with 3 column volumes of Elution Buffer followed by equilibration with at least 3 column volumes of Binding/Wash buffer. Columns can be regenerated up to 10 times without significant loss of binding capacity.





### 3.2 Cleaning-in-place

In general, **rProtein A Beads 4FF** are well suited for reuse a number of times. When precipitation and protein aggregation cause the loss of velocity and combined loads, you need to clean the medium.

#### Remove the precipitation or denatured protein

Wash the column with 2 column volumes 6M guanidine hydrochloride solution. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

#### Remove the hydrophobically bound protein

Wash the column with 3-4 column volumes 70% ethanol or 2 column volumes 0.1% non-ionic detergent. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

## 4. Troubleshooting

| Problem                                                                                            | Possible Cause                                               | Solution                                                                                                             |
|----------------------------------------------------------------------------------------------------|--------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| The flow rate of the column is very low.                                                           | The sieve plate is blocked.                                  | Clean or replace the sieve plate.                                                                                    |
|                                                                                                    | Column is clogged.                                           | Cleaning in place(part 3).                                                                                           |
|                                                                                                    |                                                              | Filtering the sample solution by passing them through a 0.22μm or 0.45μm filter.                                     |
| The curve is not stable during sample purification                                                 | Tiny air bubbles from buffer or sample.                      | De-gas buffers and samples. Do not allow the column to dry.                                                          |
| A considerable amount of sample has been loaded, but no specific antibody of interest is detected. | The concentration of antibody of interest is very low.       | Purify the antibody using the specific antigen coupled to a beads (i.e., PabPur Sulfalink Beads, Cat. No. SA018005). |
|                                                                                                    | The antibody is sensitive to low-pH elution buffer           | Neutralize the eluted fractions with Neutralization Buffer immediately after elution.                                |
|                                                                                                    | The IgG subclass does not bind to protein A.                 | Try other affinity chromatography media to purify the antibody, such as rProtein G Beads or rProtein A/G Beads.      |
| The recovery rate gradually decreases.                                                             | The sample is overloaded.                                    | Reduce the loading volume.                                                                                           |
|                                                                                                    | The column is too dirty and the binding capacity is reduced. | Cleaning in place(part 3).                                                                                           |

## 5. Related Products

| Product        | Cat. No. | Size   |
|----------------|----------|--------|
| HiPur A 4FF    | SA015C20 | 20 ml  |
| HiSelect A 4FF | SA015C47 | 4.7 ml |

