



rProtein A Beads

Index

1. Product Description.....	1
2. Purification Procedure.....	2
3. Maintenance.....	2
4. Troubleshooting.....	3
5. Related Products.....	3

1. Product Description

Smart-Lifesciences **rProtein A Beads** is an affinity chromatography medium designed for the purification of monoclonal antibody, polyclonal antibody and Fc tag fusion protein. The recombinant protein A is coupled to 4% sepharose gel. This coupling method can improve the binding capacity for immunoglobulin. The static binding capacity of **rProtein A Beads** is greater than 40 mg human IgG/ml settled beads. The dynamic binding capacity will vary depending on several factors such as target antibody, flow rate etc. The characteristics of **rProtein A Beads** are summarized in Table 1 .

Protein A, a bacterial cell wall protein isolated from *Staphylococcus aureus*, combines mammalian IgGs mainly through Fc regions. Recombinant protein A contains five high affinity IgG binding domains with other non-essential domains removed to reduce nonspecific binding.

Table 1. Characteristics of **rProtein A Beads**

Item	Description
Matrix Spherical	4% agarose
Ligand	recombinant protein A
Static Binding Capacity	>40mg human IgG/ml medium
Particle size	45-165um
Maximum Pressure	0.1MPa, 1bar
pH	3-10
Storage Solution	1×PBS containing 20% ethanol
Storage Temperature	2℃-8℃

Table 2. Relative binding strengths of antibodies from various species to protein G and protein A as measured in a competitive ELISA test.

Species	Subclass	Protein A binding	Protein G binding
Human	IgA	variable	—
	IgD	—	—
	IgE	—	—
	IgG1	++++	++++
	IgG2	++++	++++
	IgG3	—	++++
	IgG4	++++	++++
	IgGM	variable	—
Avian egg yolk	IgY	—	—
Cow		++	++++
Dog		++	+
Goat		—	++
Guinea pig	IgG1	++++	++
	IgG2	++++	++
Hamster		+	++
Horse		++	++++
Koala		—	+
Llama		—	+
Monkey(rhesus)		++++	++++





Mouse	IgG1	+	++++
	IgG2a	++++	++++
	IgG2b	+++	+++
	IgG3	++	+++
	IgM	variable	—
Pig		+++	+++
Rabbit	no distinction	++++	+++
Rat	IgG1	—	+
	IgG2a	—	++++
	IgG2b	—	++
	IgG3	+	++
Sheep		+/-	++

++++=strong binding; ++= medium binding; —=weak binding or no binding

2. Purification Procedure

2.1 Buffer Preparation

Water and chemicals used for buffer preparation should be of high purity. It is recommended filtering the buffers by passing them through a 0.22 or 0.45 µm filter before use.

Binding/Wash Buffer: 0.15 M NaCl, 20 mM Na₂HPO₄, pH 7.0

Elution Buffer: 0.1 M glycine, pH 3.0

Neutralization Buffer: 1 M Tris-HCl, pH 8.5

2.2 Sample Preparation

To insure that proper ionic strength and pH are maintained for optimal binding, it is necessary to dilute serum samples, ascite fluid or cell culture supernatant at least 1:1 with Binding/Wash Buffer. Alternatively, the sample may be dialyzed overnight against Binding/Wash Buffer.

2.3 Column Purification

- 1) Mix the slurry by gently inverting the bottle several times to completely suspend the **rProtein A Beads**. Close the column outlet leaving the net covered with packing buffer. Transfer the slurry to the column.
- 2) Allow the resin to settle down and the buffer to drain from the column. Add 5 column volumes Binding Buffer onto the column to equilibrate the beads.
- 3) Apply the sample to the column and drain the flow-through. Collect the flow-through for measuring the binding efficiency to the beads, i.e. by SDS-PAGE.
- 4) Wash the column with 10-15 column volumes Wash Buffer or until the absorbance of the effluent at 280 nm is stable.
- 5) Elute the sample with 5-10 column volumes Elution Buffer. Collect the eluate containing the target immunoglobulin and immediately neutralize to pH 7.4 with Neutralization Buffer (1/10 volume of total eluate).
- 6) Equilibrate the column with 3 column volumes of Binding Buffer , 5 column volumes of distilled water and 5 column volumes of 1XPBS containing 20% ethanol. Finally store the resin with 1XPBS containing 20% ethanol at 4°C.

2.4 Analysis

Identify the fractions using UV absorbance, SDS-PAGE, or western blot.

3. Maintenance

3.1 Regeneration

Regenerate the column by washing the resin with 3 column volumes of Elution Buffer followed by equilibration with at least 3 column volumes of Binding/Wash buffer. Columns can be regenerated up to 10 times without significant loss of binding capacity.

3.2 Cleaning-in-place

In general, **rProtein A Beads** are well suited for reuse a number of times. When precipitation and protein aggregation cause the loss of velocity and combined loads, you need to clean the medium.





Remove the precipitation or denatured protein

Wash the column with 2 column volumes 6M guanidine hydrochloride solution. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

Remove the hydrophobically bound protein

Wash the column with 3-4 column volumes 70% ethanol or 2 column volumes 0.1% non-ionic detergent. Finally wash the column with 5 column volumes 1XPBS (pH7.4).

4. Troubleshooting

Problem	Possible Cause	Solution
The flow rate of the column is very low.	The sieve plate is blocked.	Clean or replace the sieve plate.
	Column is clogged.	Cleaning in place(part 3).
		Filtering the sample solution by passing them through a 0.22µm or 0.45µm filter.
The curve is not stable during sample purification	Tiny air bubbles from buffer or sample.	De-gas buffers and samples. Do not allow the column to dry.
A considerable amount of sample has been loaded, but no specific antibody of interest is detected.	The concentration of antibody of interest is very low.	Purify the antibody using the specific antigen coupled to a beads (i.e., PabPur Sulfonink Beads, Cat. No. SA018005).
	The antibody is sensitive to low-pH elution buffer	Neutralize the eluted fractions with Neutralization Buffer immediately after elution.
	The IgG subclass does not bind to protein A.	Try other affinity chromatography media to purify the antibody, such as rProtein G Beads or rProtein A/G Beads.
The recovery rate gradually decreases.	The sample is overloaded.	Reduce the loading volume.
	The column is too dirty and the binding capacity is reduced.	Cleaning in place(part 3).

5. Related Products

Product	Cat. No.	Size
rProtein A Beads	SA012005	5 ml
	SA012025	25 ml
	SA012100	100 ml
	SA012500	500 ml
	SA01201L	1 L
	SA01210L	10 L
AbPur rProtein A Kit	SA012K03	3T
rProtein A Beads 4FF	SA015005	5 ml
	SA015025	25 ml
	SA015100	100 ml
	SA015500	500 ml
	SA01501L	1 L
	SA01510L	10 L
AbCap A 4FF	SA015C11	1 X 1 ml
	SA015C51	5 X 1 ml
	SA015C15	1 X 5 ml
	SA015C55	5 X 5 ml
	SA015CS	3X1 ml+1X5 ml

